REUSABLE CRANIAL PERFORATORS INSTRUCTIONS FOR USE

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Japanese/日本語

Chinese/中文

**Reusable Cranial Perforators
"R" Series**

**IMPORTANT INFORMATION
Read Before Using**

Model 200-141.....14/11mm DGR-I

Model 200-151.....14/11mm DGR-II

Model 200-152.....11/7mm DGR-II

Model 200-171.....14/11mm DGR-0

**IMPORTANT
INFORMATION
Read before using
FAILURE TO READ
AND FOLLOW THESE
INSTRUCTIONS MAY
RESULT IN SERIOUS
PATIENT INJURY,
DEATH OR PRODUCT
MALFUNCTION.**

Please note: These instructions include the use, care and cautions of the DGR-I, DGR-II and DGR-0 reusable perforators. Use caution when selecting appropriate model for use in specific skull areas and skull conditions.

CAUTION: ACRA-CUT cranial perforators are very sharp, therefore cut and respond differently than other brand perforators. To ensure safe and reliable performance, carefully read and follow instructions. Failure to follow instructions for use may result in serious patient injury, death or product malfunction.

**DGR-I and DGR-0
Perforators**

All Model DGR-I and DGR-0 reusable perforators are designed to automatically release and stop upon penetration of bone **at least 3mm thick**. Caution should be taken to avoid use on the temporal and suboccipital areas and on children and/or infants.

DGR-II Perforators

All Model DGR-II reusable perforators are for use on thin skull/skull areas such as pediatric, temporal, and suboccipital areas. They are designed to automatically release and stop upon penetration of bone **as thin as 1mm**. Caution should be taken to avoid excessive force when drilling thin skulls to prevent pushing the entire unit through the skull.

INSTRUCTIONS

1.0 Remove Drill from case:

With the ACRA-CUT letters facing you, place your left index finger and thumb above and below the arrow on the locking cover. Hold the case with your right hand, and with your left hand pull the cover to the left and lift. (See Picture 1)

2.0 Disassemble Unit:

- 2.1 Align the small external hole on the drill body with the hole on the inside of the drill by rotating the shank. Once the 2 holes are aligned, you will be able to see through the drill.
- 2.2 Place the drill body in the sterilizing case with the aligned holes over the pin and the cutting edges extending off the case. (See Picture 2) Close the locking cover to secure drill (reverse instructions 1.0).
- 2.3 Rotate the outer drill counter-clockwise until it becomes disengaged from the shank. Open the cover and remove drill body from the case.
- 2.4 Assembly is the opposite procedure.
- 2.5 The drive section is lubricated for life and must not be disassembled.

CAUTION:

All cutting tip components must be assembled using only the same serial numbered parts. In addition, parts between different size or model ACRA-CUT units **are not** interchangeable.

3.0 Assembly/Disassembly of Cutting Edges:

- 3.1 **14/11mm Reusable Perforators** (See Picture 3)
 - 3.1.1 Insert cam end of inner drill (a) into top of outer drill (b).
 - 3.1.2 Thread assembled parts (a/b) onto drive shank (c).
 - 3.1.3 Tighten per paragraph 2.0
- 3.2 **11/7mm Reusable Perforators** (See Picture 4)
 - 3.2.1 Insert cutting end of inner drill (a) into the base of outer drill (b).
 - 3.2.2 Slide large ring (c) over the cam end of (a) rotating until the key on (c) fits into the slot on (b).
 - 3.2.3 Place small ring (d) over the front of the (a/b) assembly and screw the large ring (c) and the small ring (d) together.
 - 3.2.4 Thread assembled parts (a/b/c and d) onto drive shank (e).
 - 3.2.5 Tighten per paragraph 2.0.

4.0 Disinfection/Manual Cleaning Instructions: Machine Cleaning is not recommended.

- 4.1 Prepare Enzyme Detergent per manufacturer's instructions.
 4.2 Disassemble perforator per ACRA-CUT's Instructions for Use (Per Paragraph 2.0)
 4.3 Perform cleaning per the steps listed in the following table.

Step	Solution	Time Minutes	Temperature	Instruction
1	Hospital Grade Enzymatic Detergent	14-15	Room Temperature	Immerse and soak for required time.
2	Hospital Grade Enzymatic Detergent	As Required	Room Temperature	Clean thoroughly - Scrub all external surfaces with a soft bristle brush until all visible soil has been removed. It is important to make sure that the flutes are effectively cleaned. Use a small diameter brush or pipe cleaner to clean cannulation holes. Inspect for visible soil on exposed surfaces.
3	Water	2-3	Warm, as delivered from hot water tap	Rinse thoroughly for required time immediately after step 2.
4	Air	As Required	Ambient	Allow to air dry in clean area. Blow lumens with clean air using filtered air source or syringe.

5.0 Sterilization:

- 5.1 Perforator must be manually cleaned prior to sterilization. See Paragraph 4.0 for recommended procedure.
- 5.2 Full Cycle Steam Guidelines are as follows:

Cycle #	Cycle Type	Temp (C°)	Exposure Time (Minutes)	Dry Time (Minutes)
1	Pre-Vac, wrapped	134	4	20
2	Gravity, wrapped	121	30	30

- 5.3 The use of other sterilization methods is outside of ACRA-CUT's responsibility. If the described equipment, materials and chemicals are not available, it is the responsibility of the user to validate their method accordingly.
- 5.4 Frequent reprocessing has few effects on these instruments. The end of the instrument's working life is normally determined by wear and damage through use.

6.0 Preoperative Tests:

Perform the following test before each perforation.

Protect gloved fingers with gauze while performing tests. If either test does not perform as indicated, take special action to prevent further use of the perforator and immediately return to ACRA-CUT, Inc. To test the mechanical integrity of the assembled unit:

- 6.1 Hold the drive shank securely while pressing the fluted components (inner and outer drill) toward each other. A smooth, positive spring action should be felt.
- 6.2 While holding the drive shank, rotate the fluted portion (inner and outer drill). These parts should turn smoothly.

7.0 General Operating Procedures:

- 7.1 The drive shank fits a standard Hudson coupler.
- 7.2 Power unit must be used only in the CLOCKWISE direction, or forward mode.
- 7.3 Perforator can be used at low speeds and in hand-operated drilling devices.

7.4 Recommended speed for power drilling is up to 800 RPM. Power unit (standard 1000 RPM) should be first set at 60 PSI and increased as needed to no more than 80 PSI.

7.5 To start perforation, place the nonskid tip against the skull, compress and rotate drill by hand until drive mechanism is seated in the driver. Then drill at chosen speed.

7.6 Perforator will stop when penetration is made or when surgeon releases pressure on the unit.

7.7 Shelf Formation:

7.7.1 The shelf formed by DGR-I and DGR-0 perforators will be well defined. The DGR series of perforators also form a ridge of bone at the bottom of the through hole. (See Figure 1)

7.7.2 The shelf formed by DGR-II perforators slants more than the shelf formed by DGR-I and DGR-0 perforators. This feature allows for safe use on skull/skull areas as thin as 1mm. The DGR series of perforators also form a ridge of bone at the bottom of the through hole. (See Figure 2)

- 7.8 The Dura-Guard “R” Series tip forms a pad of bone at the bottom of the burr hole. This pad protects the dura from the cutting edges and pushes the dura away from the skull. After irrigation and suctioning, the pad is easily removed with a pick or elevator type instrument, such as the ACRA-CUT Model 800-020. If it is suspected that the dura is adhering to the bone, EXTRA CAUTION must be taken or dura may be cut. (See Figure 3)
- 8.0 WARNING:**
Failure to heed the warnings or to follow the procedures as described could result in serious patient injury, death or product malfunction.
- 8.1 Surgeon should always be aware of the possibility that the perforator may not disengage.
- 8.2 Excessive hand pressure must be avoided to preclude penetration of the dura and/or brain. Allow the perforator to do the cutting.
- 8.3 Reduce pushing force and observe caution as the perforator nears the point of penetration as nicking, cutting, or tearing the dura can occur.
- 8.4 Hold power unit/perforator perpendicular (90 degrees) to the skull at the point of penetration and do not “rock” to speed cutting action. Failure to observe may cause product malfunction, patient injury or dural lacerations.
- 8.5 Extreme caution must be taken when perforating in areas of the skull that have 1mm or greater variations in bone thickness. Drill may nick or cut dura or brain (similar effect as not being 90 degrees).
- 8.6 If conditions such as adherent dura, high intracranial pressure or other underlying abnormalities are present in the area of the penetration, the perforator may cut the dura or brain.
- 8.7 Speed for power drilling should not exceed 800 RPM. See Paragraph 7.4 for procedure.
- 8.8 Although the perforators are designed to release safely and reliably (DGR-I/DGR-0 units in skull/skull areas at least 3mm thick; DGR-II units in skull/skull areas as thin as 1mm) caution must be observed when perforating thin skull areas, such as temporal bone, infants, children, elderly, or diseased bone, since skull consistencies and thicknesses can vary and the dura could be cut.

- 8.9 Do not use perforator in or near previously drilled areas.
- 8.10 Perforator must be disassembled, cleaned and sterilized between each use. See Paragraphs 2.0, 3.0, 4.0 and 5.0. Failure to follow these instructions could cause a product malfunction and subsequent serious patient injury or death.
- 8.11 Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress corrosion cracking.
- 8.12 The user is responsible for the actual processing achieving the desired results with the equipment, materials and staff employed in the processing facility. Usually, this requires validation and routine monitoring of the method.
- 8.13 Machine cleaning of this device is not recommended.
- 9.0 Resharpening:**
- 9.1 The manufacturer is the only authorized party to resharpen or service the unit. ACRA-CUT, Inc. assumes no liability if the sealed unit (shank end) has been disassembled or if the resharpening or servicing is performed by any party other than ACRA-CUT, Inc.
- 9.2 Quarterly no-charge inspections and recertification of the perforator at the ACRA-CUT, Inc. factory are strongly recommended.
- 9.3 For resharpening, or for a no-charge inspection, return perforator to:
ACRA-CUT, Inc.
989 Main Street
Acton, Massachusetts 01720
USA
- 9.4 Unit will be returned by the same method as shipped. Inspection/resharpening will take no more than (3) working days.

CAUTION: Federal (USA) law restricts this device to use by or at the order of a physician.

ACRA-CUT products are manufactured under one or more of the following patents:

U.S. Patent No. 4,600,006;

U.S. Patent No. 4,699,550;

U.S. Patent No. 4,803,982;

U.S. Patent No. 4,884,571;

U.S. Patent No. 4,951,690;

U.S. Patent No. 5,007,911;

U.S. Patent No. 5,135,532;

U.K. Patent No. 2,154,484;

U.K. Patent No. 2,181,076;

U.K. Patent No. 2,184,374;

U.K. Patent No. 2,208,362;

DE Patent No. 3,503,098;

DE Patent No. 3,624,860;

DE Patent No. 3,890,886;

JP Patent No. 2,607,658.


Other U.S. and foreign patents pending.

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ACRA-CUT Inc.
MT Promedt Consulting GmbH